

# OPTO CLEAN

**CLEANING CONCENTRATE FOR OPTICS**  
*universal for all spectacle frames and glasses*

## Description

*OPTO CLEAN* is an approved concentrate for the cleaning of glasses and frames in ultrasonic baths. Also suitable for the cleaning of:

- coated glasses made from plastic materials,
- framed and other precision optical pieces.

The advantages of *OPTO CLEAN* are:

- removes gently grease and sweat deposits, fingerprints,
- no damages for glasses and frames also after multiple application,
- brightening for most metallic frames,
- highly economical due to the concentrated formula.

## Application and dosage

Dilution ratio: dilute with water at about 1:20 – 1:50 (2 - 5 vol%). Depending on the degree of contamination, clean in ultrasonic bath for approx. 1 – 10 min. Rinse the cleaned pieces in running water and dry.

## Safety recommendations

*OPTO CLEAN* is not classified as hazardous according to Regulation (EC) No 1272/2008 [GHS]. Observe the hints indicated in the Safety Data Sheets. Always handle chemicals with care.

## Physical-chemical characterisation

Density: 1.02 g/ml, pH 9.5 – 10 in state of delivery.

Ingredients according to Annex VII, A, EC-Regulation 648/2004 (detergents): <5% anionic surfactants, <5% non-ionic surfactants, perfumes.

## Disposal

The surfactants in our product meet the criteria for biodegradation as laid down in Annex III of the Regulation (EC) No 648/2004 on detergents.

The cleaning bath can be fed into the sewage system after use; check the contamination brought into the used solution. Disposal of the concentrate: eventual neutralisation (observe local sewage regulations).

European waste code: 20 01 30, „detergents other than those mentioned in 20 01 29“.

## Volumes, storage and transport

Available volumes: 1 litre PE-bottle; 10 litre and 25 litre HDPE-can.

Store in closed original container at a temperature between +5°C and +30°C, protected from heat and direct solar radiation.

Shelf life: 5 years from date of production (see stamp on label).

Classification for all means of transport: no hazardous material.



**Elma Schmidbauer GmbH**

Gottlieb-Daimler-Str. 17, D-78224 Singen

e-mail: [info@elma-ultrasonic.com](mailto:info@elma-ultrasonic.com)

12.03.2018 / (GB) Version 2.0

Tel.: +49 7731 882-0

Fax: +49 7731 882-266

[www.elma-ultrasonic.com](http://www.elma-ultrasonic.com)